EN FR DE IT ES TR KR

# BONDFILL SB

## EN

### INSTRUCTIONS

**IMPORTANT: READ ALL INSTRUCTIONS THOROUGHLY BEFORE USE. KEEP THIS LEAFLET AND REFER TO IT PERIODICALLY.**

**Dental Adhesive Resin for Filling Polymer-based Pit and Fissure Sealant, Class 1**

### Bondfill SB

**FOR DENTIST USE ONLY**

**1. What is Bondfill SB?**  
Bondfill SB is a self-cure, self-adhesive resin for dental restorations. Bondfill SB is chemically activated with extra-oral. It contains 4-META/MMA-TBB. Though it requires no bonding agent, Bondfill SB forms bonds to enamel, dentin, and restorative materials without gapping and discoloration. It is also designed to absorb insect, complicated complex external stresses and resist wear. The Powder contains PMMA and TMPMT reactive organic filler which creates a tough restoration that is hard yet still resilient. The inorganic filler loading is 1.2 vol% (as powder) with a particle size range of approx. 0.02µm.

**Indications:**  
Bondfill SB is indicated wherever resilience is required, such as:  
a. Restoration of tooth wear, such as attrition, abrasion, abrasion, wedge-shaped defects.  
b. Secondary caries at a restorative margin.  
c. Repair of fractured prostheses and damaged restorations.  
d. Sealing incipient pit and fissure caries.

**Note:** Bondfill SB is not indicated to use for the restoration of occlusal surface.

### 2. Contents of the kit

#### 2.1 Bondfill SB Kit:

- ①Base ----- 1 (8mL)
- ②Catalyst V ----- 1 (0.7mL)
- ③Powder Light (approx. Shade A2) ----- 1 (3g)
- ④Powder Medium (approx. Shade A3.5) ----- 1 (3g)
- ⑤Teeth Primer ----- 1 (3mL)
- ⑥Sponge (L-S) ----- 1
- ⑦Dispensing Stand ----- 1
- ⑧Dispensing Cups ----- 1
- ⑨Brush Handle (Bent) ----- 1
- ⑩Brush Tips (Brush-dip L) Label color: Pink ----- 1
- ⑪Brush Tips (Brush-dip LL) Label color: Purple ----- 10

**Attached documents :** Instructions / Pictorial Instruction Card

- 2.2 Refills:**
- Base ----- 8mL
  - Catalyst V ----- 0.7mL
  - Powder Light (approx. Shade A2) ----- 3g
  - Powder Medium (approx. Shade A3.5) ----- 3g
  - Powder Ceramic (approx. Shade A4) ----- 3g
  - Powder Opaque (approx. Shade A3 Opaque) ----- 3g
  - Teeth Primer ----- 3mL
  - Sponge (L-S) ----- 1
  - Dispensing Stand ----- 1
  - Dispensing Cups ----- 40
  - Brush Handle (Straight) ----- 1
  - Brush Handle (Bent) ----- 1
  - Brush Tips (Brush-dip S) Label color: Green ----- 1
  - Brush Tips (Brush-dip L) Label color: Pink ----- 10
  - Brush Tips (Brush-dip LL) Label color: Purple ----- 10

**Additional documents :** Instructions / Pictorial Instruction Card

①Base -----	1 (8mL)
②Catalyst V -----	1 (0.7mL)
③Powder Light (approx. Shade A2) -----	1 (3g)
④Powder Medium (approx. Shade A3.5) -----	1 (3g)
⑤Teeth Primer -----	1 (3mL)
⑥Sponge (L-S) -----	1
⑦Dispensing Stand -----	1
⑧Dispensing Cups -----	1
⑨Brush Handle (Straight) -----	1
⑩Brush Tips (Brush-dip L) Label color: Pink -----	1
⑪Brush Tips (Brush-dip LL) Label color: Purple -----	10

**3. Compositions.**  
Base: Methacrylate monomers, others.  
Powder: Poly(methyl methacrylate, others.  
Catalyst V: Partially oxidized tri-n-butylborane, others.  
Teeth Primer: Methacrylate monomers, Aceton, Water, others.

#### 4. Precautions :

**4.1 Safety**  
This product contains substances that may cause allergic reactions. Read the following points below.

- ①**Sensitivity**  
Bondfill SB should not be used by clinicians or on patients who are sensitive to methacrylates or acetone. If signs of irritation such as redness appear, stop using Bondfill SB immediately and consult a physician.  
②**Avoid contact**  
Avoid contacting soft tissues, skin or eyes by Bondfill SB components. The use of rubber dam is recommended. Contaminated skin or mucosa should be wiped immediately with alcohol and then rinsed thoroughly with copious amounts of running water. If Bondfill SB components enter the eye, rinse immediately and thoroughly with running water and have the patient examined by an ophthalmologist. Clinicians should use dental gloves.  
③**Avoid ingestion**  
Avoid accidental ingestion while applying and rinsing this material. If ingested, the patient must seek medical attention.  
④**Flammability**  
Catalyst V, Base, and Teeth Primer are highly volatile and flammable. Do not store where they may be exposed to open flame.  
⑤**Avoid inhalation**  
Handle Bondfill SB in a well ventilated area (air to be replaced several times per hour). In case of inhalation of the vapor, move for fresh air.

**6 Applications**  
Use Bondfill SB only for the applications recommended in these instructions. If necessary, for deep cavities with only a minimal layer of dentin, the pulp must be indirectly capped liquid. Don't use in conjunction with other materials. Brush Tips is for one-time use only. Do not reuse and should be discarded after using.

#### 4.2 Precautions on Catalyst V

**Catalyst V reacts with air and water to generate heat and lose activity. Please abide by the following.**

- ①**Storage conditions**  
Avoid high temperature, high humidity and direct sunlight. The repeated temperature changes may cause the syringe to aspirate air, thereby shortening the Catalyst V's shelf-life. \*After prolonged storage, the first drop of the Catalyst V may be inactive due to exposure to air, though the rest of the material remains active. Express and dispose of a drop before making the mix.  
\*The syringe is made of glass, therefore it must be handled with care to prevent shock, dropping, and other physical damage.  
②**Cap closure**  
The cap simply slides on and off. Recap the syringe immediately after each use. Air (oxygen and humidity) deactivates the Catalyst V. Do not leave the cap off during the bonding procedure.  
③**Expressing the Catalyst V**  
Bondfill SB does not come out of the syringe because the tip is blocked or the screw is tight, do not try too hard to turn it. The contents may splash as the syringe breaks.  
④**After use**  
Unscrew the syringe screw two turns counter-clockwise after each use to relieve pressure on the Catalyst V. (Pressure buildup can cause leakage of the Catalyst V or even crack the syringe.)  
⑤**Clean spilled Catalyst V immediately with wet towels**  
The Catalyst V reacts with oxygen. If absorbed by a flammable material, it may raise the temperature enough to cause smoldering. If the Catalyst V is spilled, wipe it up immediately with a WET (not dry) disposable towel. Then rinse the towel to deactivate the Catalyst V thoroughly in running water.  
⑥**Cleaning syringe**  
Wipe the tip of the syringe with a dry gauze after each use to prevent residue buildup. Then rinse the gauze with water to kill any remaining activity. Buildup of the residue may prevent the cap from seating properly.

**4.3 Storage**  
**Please take the following precautions to maintain the quality.**

①**Storage conditions**  
Store Catalyst V, Base and Teeth Primer in a cool (1°C - 30°C) and dark location. High temperature, high humidity and direct sunlight will shorten their shelf-life.  
②**Volatility**  
Teeth Primer and Base are highly volatile. Recap the bottle immediately.  
③**Contamination**  
To avoid contamination, do not mix the components with other materials. Firmly recap the containers immediately after use. Dispose of any excess of mixed materials after use. Once dispensed do not return materials to the containers. Both Dispensing Cups and Brush Tips are for one-time use. Dispose of them after use.  
④**Shelf life**  
3 years in a room temperature.  
\*See outer package for expiration date. Do not use after expiration date.

**4.4 Disposal**  
①Soak the empty syringe of Catalyst V in water overnight. Dispose it appropriately according to the law about the waste.  
②Dispose of empty container or package except Catalyst V according to the law about the waste.

**As in any dental treatment, the patient's individual constitution and the unique requirements of clinical case at hand must be considered before selecting materials and conditions for use.**

#### 5. How to use Bondfill SB

**a. Surface preparation**  
① **Cleaning**  
Remove contaminants and stains using a polishing brush or oil-free, fluoride-free pumice to expose fresh tooth surface. Rinse thoroughly and dry. A rubber dam is highly recommended to reduce the chance of contamination by saliva, breath humidity or blood.  
② **Priming**  
Dispense Teeth Primer into one of the small recesses on the Dispensing Stand. Using a pledget, keep the surface to be bonded saturated with a copious amount of Teeth Primer for approximately 20 seconds. Care should be taken not to contact the gingiva with the pledget. Then air dry. Rinsing is not needed.

**Standard application time**

Surface	Teeth Primer	Enamel Etchant Gel	Dentin Etchant Gel
Enamel	Approx. 20 sec	30 sec	5-10 sec
Dentin	-----	-----	5-10 sec
After application	Dry without rinsing	Rinse and dry	-----

**3) Metal Surface**  
Rinse the metal surface thoroughly and dry. Applying the Enamel Etchant Gel helps remove surface contaminants. Air-abrade the metal surface with 50-micron aluminum oxide. If the surface is precious metal, apply a metal primer such as V-PRIMER on the metal surface to be bonded in order to increase bond strength and durability. V-PRIMER is sold separately by Sun Medical Co., Ltd. For application details refer to the instructions for the metal primer.

**4) Ceramic (porcelain and zirconia) surface**  
Depending on the case at hand, clean the porcelain/zirconia surface with abrasive wheels and points or a prophy cup with fluoride-free, oil-free pumice. Rinse thoroughly with water and dry. Applying the Enamel Etchant Gel helps remove surface contaminants. An ultrasonic cleaner is also effective. Then, apply a porcelain primer such as Porcelain Liner M / Super-Bond Universal Ceramic Primer on the porcelain/zirconia surface to be bonded in order to increase bond strength and durability. Porcelain Liner M / Super-Bond Universal Ceramic Primer is sold separately by Sun Medical Co., Ltd. For application details refer to the instructions for the porcelain primer.

**b. Preparation of Bondfill SB**  
Place two Dispensing Cups into two of the three holes in the Dispensing Stand. Dispense 2-3 drops of the Base into one of the Dispensing Cups. Hold the Catalyst V syringe vertically, and turn the screw to dispense a drop into the Base. Stir lightly with a brush. This mixture is called "activated liquid." The activated liquid gradually deactivates. Use it within 3 minutes after preparation.

#### Activated Liquid

Intended Use	Base	Catalyst V	Operation time
Restoration and Build-up	2-3 drops	1 drop	3 minutes after mixing
Pit and Fissure Sealant	2 drops	-----	-----

#### 2. Restoration/Build-up/Pit and Fissure Sealant

- Dispense the appropriate amount of Powder into the other Dispensing Cup.
- Wet the surface to be bonded with the activated liquid.
- Soak the brush with ample activated liquid, then dip the brush into the Powder to form a ball of powder on the tip. Brush the powder ball onto the pre-wet surface to be bonded.
- Repeat the procedure until the restoration or the build-up is completed. When you repeat the procedure, clean the brush tip with gauze before you dip it again.

**d. Setting and finishing**  
Wait approximately 10 minutes before finish and polish. If complete setting is desirable, wait 10 points until the next day.  
\*Use carbide burs or super-fine diamond points for finishing.  
\*Use white silicone points for polishing.

**Working time and Setting time of Sealant**

Item	Specification	Standard
Working time	≤5 minutes	ISO6874(2005)
Setting time	≥40 seconds	ISO6874(2005)

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