# Metafil Bulk Fill

EN

## INSTRUCTIONS

**IMPORTANT:**  
 READ ALL INSTRUCTIONS THOROUGHLY BEFORE USE.  
 KEEP THIS LEAFLET AND REFER TO IT PERIODICALLY.

Resin-based Dental Restorative Material

## Metafil Bulk Fill

Type1, Class2, Group1

**FOR DENTIST USE ONLY**

### 1. Product Overview

Metafil Bulk Fill is a resin-based light-cure, radiopaque, flowable restorative material. The content of inorganic fillers is 47.6 vol % (70.0 wt %). The particle size of the inorganic fillers is between 0.01-8.00  $\mu\text{m}$ .

#### Intended use:

- Metafil Bulk Fill is intended to be used as a resin-based dental restorative material.

#### Indication:

- Direct/indirect posterior restoration including occlusal surfaces.

#### Contraindications:

Patients who have a history of severe allergic reactions to this product, methacrylate monomers, or acrylate monomers.

#### Side effects:

In rare cases, components of Metafil Bulk Fill may lead to sensitization. In such cases, Metafil Bulk Fill should no longer be used.

### 2. Components

Metafil Bulk Fill ----- 2.5mL (4.8g)

<Low-Flow>

Shade: Universal, A1, A2, A3

<Medium-Flow>

Shade: Universal, A1, A2, A3, Blue

<Needle Tip in Low-Flow package>

18G Needle Tip (with green ring) ----- 5pcs

<Needle Tip in Medium-Flow package>

19G Needle Tip (with pink ring) ----- 5pcs

■ Attached documents: Instructions, Shade Guide

### 3. Compositions

Metafil Bulk Fill is a paste consisted of the following compositions.

Methacrylate monomers, Acrylate monomers, Inorganic fillers, Aromatic amine, and the others.

### 4. Precautions

This product contains substances that may cause allergic reactions. Read the following points below.

#### 1) Safety Precautions

- ① Metafil Bulk Fill should not be used by clinicians or on patients who are methacrylate monomer-sensitive. If signs of irritation such as redness appear, stop using it immediately and consult a physician.
- ② Avoid contact of Metafil Bulk Fill with soft tissues, skin or eyes. Contaminated skin or mucosa with Metafil Bulk Fill should be wiped off immediately with alcohol and then rinsed thoroughly with running water.
- ③ If Metafil Bulk Fill enters an eye, immediately rinse the eye thoroughly with running water. The patients should then be examined by an ophthalmologist.
- ④ Avoid accidental ingestion while applying this product. If ingested, the patient must seek medical attention.

#### 2) Common Precautions

- ① Dentist should wear rubber or PVC dental gloves and eyewear.
- ② Proper use of eye protection is recommended when using light irradiation devices.
- ③ Use a dental rubber dam to isolate the treatment area from contamination.
- ④ Do not mix the components with other materials. Once dispensed, do not return this product to the containers. Use Metafil Bulk Fill only for the applications recommended in this leaflet.

#### 3) Handling Precautions

- ① In the case of pulp exposure or deep cavities with close proximity to pulp, cap the pulp with calcium hydroxide etc. Do not use eugenol-containing material.
- ② Needle Tip is for one-time use only. Dispose it after each patient to avoid cross-contamination. After use, discard any excess.
- ③ Needle Tip is not autoclavable or chemically sterilizable.
- ④ Tightly close syringes with syringe cap immediately after use. Do not allow Needle Tip to remain as a syringe cap.
- ⑤ Use the Needle Tip only for this product. Use 18G Needle Tip (with green ring) only for Low-Flow and 19G Needle Tip (with pink ring) only for Medium-Flow.

### 5. Storage

#### 1) Storage conditions

- ① Store Metafil Bulk Fill at a room temperature between 1 – 30°C. Avoid high temperatures above 30°C, high humidity and direct sunlight.
- ② Shelf life: 3 years in a room temperature.
- ③ Do not use this product after the expiration date. See outer package for expiration date.

#### 2) Product lifetime

Product lifetime can be shortened depending on the circumstances of use. Read instructions of relevant products carefully before use.

### 6. Disposal

Dispose of empty containers, packaging, etc. with the normal household waste.

**As in any dental treatments, the patient's individual constitution and the unique requirements of each clinical case at hand must be considered before selecting materials and conditions for use.**

### 7. How to use

#### <Metafil Bulk Fill syringe>

- ① Remove the syringe cap by twisting it counter-clockwise. Save the syringe cap.
- ② Install the Needle Tip into the syringe certainly and make sure the Needle Tip and the syringe are connected tightly.

#### <Direct restoration>

- ① Prepare the cavity by removing caries and any existing restorations as indicated.
- ② Protect pulp by covering the dentin close to the pulp with adequate pulp capping materials.
- ③ Adhesive Bonding System:  
Apply light-cured bonding system according to the manufacturer's instructions. We recommend using bonding system manufactured by Sun Medical Co., Ltd.
- ④ Dispense Metafil Bulk Fill directly into the cavity by slowly squeezing the plunger. Cavities can be filled up to 4mm in one increment.
- ⑤ Light cure Metafil Bulk Fill. Refer to Table 1 for the curing time.  
**Note:** For deeper preparations, add more Metafil Bulk Fill by increments.
- ⑥ Dispose Needle Tip properly after use and seal the container immediately with the syringe cap. Do not store the syringe with Needle Tip in place. Store it only with the syringe cap.
- ⑦ Check occlusion with a thin articulating paper. Examine centric and lateral excursion contacts. Carefully adjust occlusion by removing material with a fine polishing diamond or stone.

Table 1 Curing time and curing depth of Metafil Bulk Fill

Type	Light intensity (mW/cm <sup>2</sup> )	Curing depth (mm)	Curing time (sec)
LED	> 1,000	4.0	10

- ⑧ Metafil Bulk Fill is polished quickly and easily to a high luster. Most resin discs or polishing disc, paste and polishing systems are acceptable.  
When Metafil Bulk Fill is used as a base/liner, please follow the directions below.

#### <Base/liner under direct restorations>

- ① Prepare the tooth in the same way as for direct restoration, step ① – ⑥.
- ② Place a surface layer of a methacrylate-based universal/posterior composite resin restorative. When Metafil Bulk Fill is used as a base/liner, surface layer composite resin material should have at least 2 mm thickness.
- ③ For working and finishing directions, follow the manufacturer's instruction for use.

#### <Indirect restoration with prosthesis>

- ① Prepare the cavity by following the step ① – ⑥ for direct restoration to place Metafil Bulk Fill as base/liner. Finish the preparation.
- ② Take an impression with agar/alginate or silicone materials.
- ③ While waiting for the fabrication of final prosthesis, properly protect the tooth surface with temporary restorative material or temporary crown. When using resin-based temporary restorative material, apply dental separating agent to the cavity wall.
- ④ Bond the prepared prosthesis with commercially available dental cement. We recommend using dental cement manufactured by Sun Medical Co., Ltd. Prime the prosthesis if indicated by the instruction of the preferred cement. We recommend dental prosthesis priming agent manufactured by Sun Medical Co., Ltd.  
**Note:** Do not use eugenol-containing temporary restorative material.

### 8. To get the best results with Metafil Bulk Fill

- ① Use 18G Needle Tip with green ring only for Low-Flow and 19G Needle Tip with pink ring only for Medium-Flow only.
- ② Allow this product to reach room temperature prior to use for proper flow if Metafil Bulk Fill is stored in a cool place.
- ③ Metafil Bulk Fill is light sensitive material. Protect from ambient or operatory lights. Reduce/turn off the light while the restorative procedure as indicated to maintain proper working time.
- ④ If the cavity is very deep, incremental restoration may be required (refer to Table 1).
- ⑤ Check the performance of the light curing device. Clean and exchange the light bulb periodically. Insufficient polymerization from the compromised irradiator may cause the lower adhesive strength.

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SUN MEDICAL CO.,LTD.