INSTRUCTIONS

IMPORTANT :
READ ALL INSTRUCTIONS THOROUGHLY BEFORE USE.
KEEP THIS LEAFLET AND REFER TO IT PERIODICALLY.

Flowable Radiopaque Hybrid Light-curing
Resin-based Dental Restorative Material

Metafil Flo

FOR DENTIST USE ONLY

Metafil Flo

1. What is Metafil Flo?

Metafil Flo is a highly-flowable, visible-light-cure, tooth-shade composite resin based on unique monomer combination and containing inorganic fillers and a unique organic filler which is reactive with the matrix. This feature assures easy handling and gives the finished restoration wear-resistance and radio-opacity. Shades are precisely matched to the Vita®-shade guide making color selection both easy and reliable. With standard polishing procedures, Metafil Flo gives outstanding life-like cosmetic results.

Typical Applications :

Metafil Flo is indicated for use mainly in smaller and complicated cavities, Class III and V cavities, and base restorations for deep Class I and II cavities, owing to its good flowability.

2. Types

Individual Syringes: Containing 1.5mL (2.6g) each
with Needle Tips 10pieces
Shades : A1, A2, A3, A3.5,
B2

Refill of the Needle Tips (50pieces/bag) is available

3. Precautions

- Metafil Flo should NOT be used either by clinicians or on patients who are methacrylate-monomer sensitive.
- Avoid contact of Metafil Flo with soft tissues, skin or eyes. Contaminated skin or mucosa with Metafil Flo should be wiped off immediately with alcohol and then rinsed thoroughly with running water. If Metafil Flo enters an eye, rinse the eye with running water immediately. The patients should then be examined by an ophthalmologist.
- Dentist should wear rubber or PVC dental gloves.
- Use Metafil Flo only for the applications recommended in this leaflet.

4. Storage

Avoid high temperatures above 30°C and direct sunlight. Refrigeration is not necessary. However, if stored in a refrigerator, then allow the material to warm to room temperature prior to use.

5. How to use Metafil Flo**Adhesive Bonding System :**

Metafil Flo is compatible with currently available dentin bonding agents including AQ Bond Plus from Sun Medical Co., Ltd. Apply the adhesive agent on the prepared surface according to the relevant instructions.

Placement and Curing :

- A. Remove the cap by turning it counter-clockwise.
 - B. Screw the Needle Tip tightly to the syringe.
 - C. Dispense Metafil Flo through the Needle Tip to the cavity previously treated with a dentin adhesive.
 - D. Light-cure for at least 30 seconds using dental visible-light curing device, remove the Needle Tip and recap the syringe.
- For best handling, use Metafil Flo at room temperature.
 - The use of an unfilled resin prior to Metafil Flo placement is unnecessary in most cases. However, should the clinician desire an unfilled resin, Metafil Flo is compatible with most bis-GMA and UDMA type unfilled resins.
 - If the cavity preparation is large and/or deep, incremental filling techniques should be employed.

Polishing :

Metafil Flo will polish quickly and easily to a high luster. Use of most discs or polishing paste systems are acceptable.

As in any dental treatment, the patient's individual constitution and the unique requirements of clinical case at hand must be considered before selecting materials and conditions for use.

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*Vita is a trademark of Vita Zahnfabrik
SUN MEDICAL CO.,LTD.