# Orthomite LC

## Light-Cure Orthodontic Adhesive

### INSTRUCTIONS FOR USE

**IMPORTANT:**  
READ ALL INSTRUCTIONS THOROUGHLY BEFORE USE.  
KEEP THIS LEAFLET AND REFER TO IT PERIODICALLY.

Light-Cure Orthodontic Adhesive

## Orthomite LC

FOR DENTIST USE ONLY

### 1. What is Orthomite LC?

Orthomite LC is a light-cure adhesive resin for orthodontic, which bonds metal, ceramic and plastic brackets to tooth surface.

Through the innovative design of the base monomer containing adhesive promoting monomer, Orthomite LC Paste can achieve strong bond to etched enamel without priming as required in most other orthodontic adhesives. Orthomite LC Paste shows excellent bond strength to metal, ceramic and plastic brackets.

### 2. Contents of the kit:

- ① Orthomite LC Paste ----- 4g x 2
- ② Orthomite LC Etchant ----- 5mL
- ③ Sponge (L•S) ----- 1
- Attached documents:
  - Instructions / Pictorial Instructions

#### Refills:

- Orthomite LC Paste ----- 4g
- Orthomite LC Etchant ----- 5mL
- Sponge (L•S) ----- 1

### 3. Precautions: Read all instructions thoroughly before use.

#### 3.1 Safety

- ① **Sensitivity**  
Orthomite LC should not be used by clinicians or on patients who are sensitive to methacrylates. If signs of irritation such as rash appear, stop using Orthomite LC immediately and consult a physician.
- ② **Avoid contact**  
Avoid contact to soft tissues, skin and eyes by Orthomite LC components. The use of rubber dam and appropriate eye protection are recommended. If contact occurs, skin or mucosa should be wiped immediately with alcohol, then rinse thoroughly with copious amounts of running water. If Orthomite LC components enter the eye, rinse immediately and thoroughly with running water and have the patient examined by an ophthalmologist. Clinicians should use dental gloves and eye protection gear.
- ③ **Be careful of acidity of Orthomite LC Etchant**  
Orthomite LC Etchant is acidic. Avoid contact with soft tissue, skin and eyes. Use appropriate isolation and evacuator to protect patients from swallowing the gel while applying and rinsing it off.
- ④ **Applications**  
Use Orthomite LC only for the applications recommended in these instructions.

#### 3.2 Storage Please take the following precautions to maintain the quality.

- ① **Storage conditions**  
High temperature, high humidity and direct sunlight will shorten the shelf-life of the components. Orthomite LC Paste should be stored IN A REFRIGERATOR (between 1- 10°C/34-50°F).  
Orthomite LC Etchant should be stored between 1- 30°C/34-86°F).

#### 3.3 To get the best bonding Results

- ① Use fluoride-free polishing paste to clean the enamel. Isolate teeth to avoid re-contamination.
- ② Do not use the brackets without mechanical retention structure in their bases.
- ③ Insufficient light output may result in incomplete cure of the adhesive paste and compromised bonding performance. Do not use the light-cure device with its output less than 300 mw/cm<sup>2</sup>.
- ④ After dispensing Orthomite LC Paste, re-cap the syringe immediately. Avoid exposure to direct light.

#### 3.4 Avoid Contamination

- ① To avoid contamination, do not mix the components with other materials. Firmly recap the syringes immediately after use.

***As in any dental treatment, the patient's individual constitution and the unique requirements of clinical case at hand must be considered before selecting materials and conditions for use.***

### 4. How to use Orthomite LC

#### ① Surface preparation

##### 1) Tooth surface

- a. Polish the surface to bond with an oil-free fluoride free pumice or paste, rinse and dry. Isolate.
- b. Using an application sponge, apply Orthomite LC Etchant to the tooth. Wait for 30 sec., rinse thoroughly with water and dry.

##### 2) Gold alloy surface

Abrade the surface by aluminum sandblast or other ways. Rinse thoroughly and dry, then brush a thin coat of precious metal primer such as V-PRIMER (Sun Medical).

##### 3) Ceramic or cured resin surfaces

Abrade the surface using aluminum sandblast or other ways. Rinse thoroughly and dry, then apply a thin coat of ceramic primer such as Porcelain Liner M (Sun Medical) and blow dry with gently air.

#### ② Apply Orthomite LC Paste

Apply a small amount of Orthomite LC Paste onto the bracket base. Recap the paste syringe immediately.

#### ③ Placement and curing

- 1) Lightly place bracket onto tooth surface.
- 2) Adjust the position and seat the bracket.
- 3) Remove excess paste around the bracket base.
- 4) Light cure the adhesive paste according to the following table.

#### Light curing conditions for different brackets (LED, Halogen)

Type of bracket	Light curing condition
Metal bracket	Both mesial and distal sides ; 10 sec. each side
Ceramic and plastic brackets	Perpendicular to bracket surface; 10 sec.

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- ①Orthomite LC Paste----- 4g x 2
- ②Orthomite LC Etchant----- 5mL
- ③Sponge (L•S) ----- 1
- Attached documents:  
Instructions / Pictorial Instructions

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**SUN MEDICAL CO.,LTD.**