EN FR DE IT ES KR

Super-Bond SEP

EN

INSTRUCTIONS

IMPORTANT : READ ALL INSTRUCTIONS THOROUGHLY BEFORE USE. KEEP THIS LEAFLET AND REFER TO IT PERIODICALLY.

Water Soluble Release Agent

Super-Bond SEP

FOR DENTIST USE ONLY

1. What is Super-Bond SEP?

Super-Bond SEP forms a thin and uniform coat and functions as an intraoral dental releaser to prevent adhesion of dental resin materials such as bonding agents, cements, to those area where bonding in undesirable. It contains hydrophillic polymers and the coated releasing film softens in water and is removed easily by washing with water.

Indications:

- When cementing a prosthesis, precoat the surfaces not to be bonded for easier cleaning.
- When preparing a temporary crown of a resin, precoat the abutment tooth for easier separation.
- In the case of temporization of a cavity protected with a resin (liner, coat, composite resin etc.), coat the resin surface before placing a resinous temporary stopping for later easier removal.

2. Contents:

Super-Bond SEP ----- 5mL

3. Precautions : Read all instructions thoroughly before use.

3.1 Safety

- ①Sensitivity** Super-Bond SEP should not be used by clinicians or on patients who has past history of sensitivity to this product. If signs of irritation such as redness appear, stop using Super-Bond SEP immediately and consult a physician.
- ②Avoid contact** Avoid contact of Super-Bond SEP with soft tissue, skin or eyes. Contaminated skin or mucosa should be wiped off immediately with alcohol and then rinsed thoroughly with copious amounts of running water. If Super-Bond SEP enters eye, immediately rinse the eye thoroughly with running water. Then the patient should be examined by an ophthalmologist. Clinicians should use dental protective gloves.
- ③Flammability** Super-Bond SEP is flammable. **KEEP AWAY FROM OPEN FLAME** in use and storage.
- ④Avoid ingestion** Do not make patients swallow Super-Bond SEP during use.
- ⑤Applications** Use Super-Bond SEP only for the applications recommended in the instructions.

3.2 To get the best results with Super-Bond SEP

- ①Storage conditions** Store Super-Bond SEP in a dark location between 15°C and 30°C. Do NOT store in a refrigerator to avoid freezing. If stored in a refrigerator and solidified, allow the material to warm to room temperature and to liquefy before use.
- ②Volatility** Super-Bond SEP is volatile. **Firmly recap the bottle immediately after use.**
- ③Avoid overrun** Apply Super-Bond SEP only to the surface to be released. If overrun onto the surface to be bonded, wash the area.
- ④Avoid Contamination** To avoid contamination, do not mix Super-Bond SEP with other materials. Once dispensed do not return materials to the containers.

As in any dental treatment, the patient's individual constitution and the unique requirements of clinical case at hand must be considered before selecting materials and conditions for use.

4. How to use Super-Bond SEP

a) Application

- Dispense a drop of Super-Bond SEP into a dish. Apply the liquid using a brush or a sponge pledget to desired surfaces. For strong releasing effect, repeat the application to form a thick coat.
- Evaporate the solvent with a gentle flow of air.
- Cap the bottle immediately after use to prevent the evaporation of the solvent.

b) Removal

- Remove resins (adhesives, temporary cements, etc.) on the Super-Bond SEP coat using a scraper, explorer, etc.
- Wash away the Super-Bond SEP until the blue color disappears. In the case it is not completely rinsed off by the water, gently rub the surface with a mechanical cleaner, such as tooth brush, moistened cotton.

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